

List of Aviv Consumable and Service Items

| Aviv Part Number | Mfr | Mfr PN | Description | Used in Aviv Model | Vendor | Telephone |
|------------------|-------|---------|---|---------------------------------|--|--------------|
| 0710-102000-8 | Osram | | XBO450W/4 | 62 and 202 | Cadmet Inc. | 610-640-1234 |
| 0710-102000-5 | Osram | 69238 | XBO150W/4 | 215 and up | Cadmet Inc. | 610-640-1234 |
| 0710-102001-5 | Osram | 69231 | XBO75W/2 | 14 (mfr 2006 and up), and 18 | Cadmet Inc. | 610-640-1234 |
| 0710-102000-9 | Ushio | 5000357 | UXL-451-0 | ATF105 | Ushio | 800-326-1960 |
| 0016-300765-0 | | | UXL-450SO MODIFIED LAMP. When ordering this lamp from Cadmet Inc., give them part number 0016-300765-0 | ATF107 | Cadmet Inc. (see note in Description column) | 610-640-1234 |
| 0710-101000-3 | Osram | 64655 | 24v 250w tungsten lamp | | Cadmet Inc. | 610-640-1234 |

For service on your Oxygen Analyzer contact: Illinois Instruments
 2401 Hiller Ridge Road
 Johnsburg, IL 60050
 Phone: 815-344-6212
 Contact: Tony DelGrosso
 Email: service2_usa@systechillinois.com
<https://www.systechillinois.com/en/support/>

For recharging your Oxygen Scrubber contact: Labclear
 Diamond Tool & Die Inc
 508 29th Avenue
 Oakland, CA 94601
 Phone: 510-534-7050 ext 202
<http://www.labclear.com/index.html>

For service on your Thermo Cube Chiller contact: Solid State Cooling
 167 Myers Corners Road
 Wappingers Falls, NY 12590
 Phone: 845-296-1300
<http://www.sscooling.com/customer-support>

For service on your Titrator contact: Hamilton Company
 Phone: 800-648-5950
<https://www.hamiltoncompany.com/support/laboratory-automation/technical-support-rno>

For Technical Service on your Aviv Spectrophotometer or Circular Dichroism Spectrometer:

Contact: Dr. Glen Ramsay
 Sundew Rising, LLC
 17 Wendy Lane
 Toms River, NJ 08753-4530
 Cell Phone: 732-278-0768
 Email: glen@sundewrising.com